




CÁMARAS DE CLIMA VARIABLE MKT

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	Intervalo de Humedad (H.R.)
 MKT 720	700	158.0x200.5x114.0	120.0x102.0x60.0	-70 a 180	-
 MKT 247	247	111.5x194.0x92.5	73.5x70.0x44.3	-70 a 180	-
 MKT 115	115	98.0x172.5x86.5	60.0x48.0x40.0	-70 a 180	-