




INCUBADORA REFRIGERADA KT

Modelo	Volumen Efectivo (L)	Dimensiones externas (bxaxp) cm	Dimensiones internas (bxaxp) cm	Rango de Temperatura (°C)	
	KT 170	163	86.0x102.5x65.5	60.0x76.5x35.5	4 a 100
	KT 115	102	86.0x71.5x65.5	60.0x45.5x35.5	4 a 100
	KT 53	53	66.0x63.5x63.0	40.0x40.0x33.4	4 a 100